

BABY-OVERHOLT LIGATURE FORCEPS

curved; serrated

13.5cm (5.4")




SKU: 040-424-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ009
KLS Martin	13-339-14
AUTHORS	BABY , OVERHOLT
Medicon	16.65.13
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ009
KLS Martin	13-339-14

AUTHORS

[BABY](#), [OVERHOLT](#)

Medicon

16.65.13

SPECIFICATIONS

[curved](#), [serrated](#)