

## BABY-OVERHOLT LIGATURE FORCEPS

curved; serrated

13.5cm (5.4")




**SKU:** 040-424-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ009
<b>KLS Martin</b>	13-339-14
<b>AUTHORS</b>	<a href="#">BABY</a> , <a href="#">OVERHOLT</a>
<b>Medicon</b>	16.65.13
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BJ009
<b>KLS Martin</b>	13-339-14

**AUTHORS**

[BABY](#), [OVERHOLT](#)

**Medicon**

16.65.13

**SPECIFICATIONS**

[curved](#), [serrated](#)