

## BAJONETT DISC RONGEUR

angled up; 3 x 9mm; 30mm

16.0cm (6.4"); CC-coated



**SKU:** 240-591-03CC

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">BAJONETT</a>
<b>SPECIFICATIONS</b>	<a href="#">angled up</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>MATERIAL</b>	<a href="#">Ceramic(CC)</a>
<b>AUTHORS</b>	<a href="#">BAJONETT</a>
<b>SPECIFICATIONS</b>	<a href="#">angled up</a>