

BAJONETT DISC RONGEUR

straight; 2 x 8mm; 30mm

16.0cm (6.4"); CC-coated



SKU: 240-590-02CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	straight