

## BAJONETT DISC RONGEUR

angled up; 2 x 8mm; 30mm

16.0cm (6.4")



**SKU:** 240-591-02

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Laminectomy Rongeurs](#)

**AUTHORS** [BAJONETT](#)

**SPECIFICATIONS** [angled up](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Laminectomy Rongeurs](#)

**AUTHORS** [BAJONETT](#)

**SPECIFICATIONS** [angled up](#)