

BAJONETT DISC RONGEUR

angled up; 3 x 9mm; 30mm

16.0cm (6.4")



SKU: 240-591-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	BAJONETT
SPECIFICATIONS	angled up



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	BAJONETT
SPECIFICATIONS	angled up