

BAJONETT DISC RONGEUR

straight; 2 x 8mm; 30mm

16.0cm (6.4")



SKU: 240-590-02

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	BAJONETT
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	BAJONETT
SPECIFICATIONS	straight