

BAJONETT PUNCH; 40° UP CUT

2mm; without ejector; thin footplate

16.0cm (6.4"); 10mm open; CC-coated




SKU: 240-090-02TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	BAJONETT
SPECIFICATIONS	40° up cut , without ejector