

BAJONETT PUNCH; 40° UP CUT

3mm; without ejector

16.0cm (6.4"); 10mm open




SKU: 240-090-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	BAJONETT
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	BAJONETT
SPECIFICATIONS	40° up cut , without ejector