

## BALL ELECTRODE; 9 FR.

90° retrograd angled

for Ø 1.9mm endoscope; 0°



**SKU:** 630-512-60

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Pediatric Resectoscopy</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Pediatric Resectoscopy</a>
<b>TYPE</b>	<a href="#">Electrodes</a>