

BECKMANN-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; 25X40mm

31.0cm (12.4")




SKU: 070-870-03

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ADSON , BECKMANN
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
AUTHORS	ADSON , BECKMANN
SPECIFICATIONS	sharp , toothed