

BECKMANN-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; 25X40mm

31.0cm (12.4")




SKU: 070-870-03

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| AUTHORS | ADSON , BECKMANN |
| SPECIFICATIONS | sharp , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Selfretaining Retractors |
| AUTHORS | ADSON , BECKMANN |
| SPECIFICATIONS | sharp , toothed |