

BECKMANN-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; 25X45mm

31.0cm (12.4")




SKU: 070-870-04

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-778-03
AUTHORS	ADSON , BECKMANN
Medicon	24.34.45
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-778-03
AUTHORS	ADSON , BECKMANN
Medicon	24.34.45

SPECIFICATIONS

[sharp, toothed](#)