

BECKMANN-ADSON SELFRETAINING RETRACTOR

blunt; 4X4 teeth; 25X25mm

31.0cm (12.4")




SKU: 070-871-01

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV271
KLS Martin	15-778-00
AUTHORS	ADSON , BECKMANN
Medicon	24.34.26
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV271
KLS Martin	15-778-00

AUTHORS

[ADSON](#), [BECKMANN](#)

Medicon

24.34.26

SPECIFICATIONS

[blunt](#), [toothed](#)