

## BECKMANN-ADSON SELFRETAINING RETRACTOR

blunt; 4X4 teeth; 25X31mm

31.0cm (12.4")




**SKU:** 070-871-02

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>AUTHORS</b>	<a href="#">ADSON</a> , <a href="#">BECKMANN</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>AUTHORS</b>	<a href="#">ADSON</a> , <a href="#">BECKMANN</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>