

## BECKMANN-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; 25X25mm

31.0cm (12.4")




**SKU:** 070-870-01

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV272
<b>KLS Martin</b>	15-776-31
<b>AUTHORS</b>	<a href="#">ADSON</a> , <a href="#">BECKMANN</a>
<b>Medicon</b>	24.34.25
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV272
<b>KLS Martin</b>	15-776-31

## AUTHORS

[ADSON](#), [BECKMANN](#)

## Medicon

24.34.25

## SPECIFICATIONS

[sharp](#), [toothed](#)