

BECKMANN-ADSON SELFRETAINING RETRACTOR

sharp; 4X4 teeth; 25X25mm

31.0cm (12.4")




SKU: 070-870-01

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV272
KLS Martin	15-776-31
AUTHORS	ADSON , BECKMANN
Medicon	24.34.25
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV272
KLS Martin	15-776-31

AUTHORS

[ADSON](#), [BECKMANN](#)

Medicon

24.34.25

SPECIFICATIONS

[sharp](#), [toothed](#)