

BECKMANN SELFRETAINING RETRACTOR

blunt; 4X5 teeth

33.0cm (13.2")


SKU: 070-691-32

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV273
KLS Martin	15-777-32
AUTHORS	BECKMANN
Medicon	24.32.06
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV273
KLS Martin	15-777-32

AUTHORS

[BECKMANN](#)

Medicon

24.32.06

SPECIFICATIONS

[blunt](#), [toothed](#)