

BECKMANN SELFRETAINING RETRACTOR

sharp; 4X4 teeth

31.0cm (12.4")




SKU: 070-690-31

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV272
KLS Martin	15-776-31
AUTHORS	BECKMANN
Medicon	24.32.04
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV272
KLS Martin	15-776-31

AUTHORS

[BECKMANN](#)

Medicon

24.32.04

SPECIFICATIONS

[sharp, toothed](#)