

BENGOLEA HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")



SKU: 040-384-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
KLS Martin	13-390-20
AUTHORS	BENGOLEA
Medicon	16.76.20
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
KLS Martin	13-390-20
AUTHORS	BENGOLEA
Medicon	16.76.20

SPECIFICATIONS

[serrated](#), [straight](#)