

## BENGOLEA HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")



**SKU:** 040-384-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-390-20
<b>AUTHORS</b>	<a href="#">BENGOLEA</a>
<b>Medicon</b>	16.76.20
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-390-20
<b>AUTHORS</b>	<a href="#">BENGOLEA</a>
<b>Medicon</b>	16.76.20

## SPECIFICATIONS

[serrated](#), [straight](#)