

BENNETT BONE LEVER

60mm; blunt


SKU: 190-580-60

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK180
KLS Martin	23-561-24
AUTHORS	BENNETT
Medicon	77.12.06



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bone Levers
Aesculap	FK180
KLS Martin	23-561-24
AUTHORS	BENNETT
Medicon	77.12.06