


BERLIN TUBING CLAMP

straight; serrated

14.5cm (5.8")



SKU: 032-006-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN908
KLS Martin	13-940-14
AUTHORS	BERLIN
Medicon	17.73.14
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN908
KLS Martin	13-940-14

AUTHORS

[BERLIN](#)

Medicon

17.73.14

SPECIFICATIONS

[serrated](#), [straight](#)