

## BEYER PUNCH

90° upwards; 2mm

shaft 9.5cm (3.8")



**SKU:** 111-502-19

**Categories:** [Open Surgery](#), [Otology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Punches</a>
<b>Aesculap</b>	OF768
<b>KLS Martin</b>	36-463-19
<b>AUTHORS</b>	<a href="#">BEYER</a>
<b>Medicon</b>	62.62.72

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Punches</a>
<b>Aesculap</b>	OF768
<b>KLS Martin</b>	36-463-19
<b>AUTHORS</b>	<a href="#">BEYER</a>

**Medicon**

62.62.72