


## BIEMER-MÜLLER APPROXIMATOR

jaw opening 5.0mm; length 7.0mm

closing force 30-40g



**SKU:** 026-540-07

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>Aesculap</b>	FD613
<b>AUTHORS</b>	<a href="#">BIEMER</a> , <a href="#">MUELLER</a>
<b>Medicon</b>	56.76.09

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>Aesculap</b>	FD613
<b>AUTHORS</b>	<a href="#">BIEMER</a> , <a href="#">MUELLER</a>
<b>Medicon</b>	56.76.09