


BIEMER-MÜLLER APPROXIMATOR

straight; jaw opening 5mm; length 9mm

closing force 30-40g



SKU: 026-540-09

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
AUTHORS	BIEMER, MUELLER
Medicon	56.74.09

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Clips, Approximators
AUTHORS	BIEMER, MUELLER
Medicon	56.74.09