


## BIEMER-MÜLLER APPROXIMATOR

angled; jaw opening 5mm; length 7.5mm

closing force 30-40g



**SKU:** 026-541-08

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>Aesculap</b>	FD614
<b>AUTHORS</b>	<a href="#">BIEMER</a> , <a href="#">MUELLER</a>
<b>Medicon</b>	56.75.09

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>Aesculap</b>	FD614
<b>AUTHORS</b>	<a href="#">BIEMER</a> , <a href="#">MUELLER</a>
<b>Medicon</b>	56.75.09