


## BIEMER-MÜLLER APPROXIMATOR

straight; jaw opening 3mm; length 6mm

closing force 30-40g



**SKU:** 026-540-06

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>AUTHORS</b>	<a href="#">BIEMER, MUELLER</a>
<b>Medicon</b>	56.74.06

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Clips, Approximators</a>
<b>AUTHORS</b>	<a href="#">BIEMER, MUELLER</a>
<b>Medicon</b>	56.74.06