

## BILIARY STENT; DOUBLE-PIGTAIL SINGLE USE

polyethylene with distal gold ring

10fr./8cm (Flaps); length: 11cm



**SKU:** PK-P108DP

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)