

## BILIARY STENT; DOUBLE-PIGTAIL SINGLE USE

polyethylene with distal gold ring

8,5fr./13cm (Flaps); length: 16cm



**SKU:** PK-P8513DP

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)