

## BIOPSY FORCEPS, FLEXIBLE

toothed spoon fenestrated, with pin

1.8mm, length: 120cm



**SKU:** PK-6BGPK1804120Y

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Biopsy Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Biopsy Forceps](#)