

## BIOPSY FORCEPS, FLEXIBLE

toothed spoon, fenestrated, with pin

1.8mm; length:: 160cm



**SKU:** PK-6BGPK1804160Y

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Biopsy Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Biopsy Forceps](#)