

## BIOPSY FORCEPS, FLEXIBLE

short-oval spoon fenestrated, with pin

2.5mm; length: 230cm



**SKU:** PK-6BOPK2504230Y

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Biopsy Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Biopsy Forceps](#)