


BIPOLAR CABLE 3.0M FLAT-CABLE PLUG

international Models; 2-pins

bipolar System; Umax: 0.8 KVp

ATTENTION: OUR BIPOLAR RESECTOSOCOPY SYSTEM CAN BE USED ONLY WITH HF-UNITS WHICH HAVE THE POSSIBILITY OF BIPOLAR CUT IN SALINE SOLUTION !!!



SKU: 287-215-30BR

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables