

BIPOLAR CABLE; FLAT-CABLE PLUG

Valleylab; international Models; 5m

unit side Ø4,0mm; Umax: 0.8 KVp

Instrument-side

Unit-side



SKU: 287-215-50

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables