

BIPOLAR CABLE; FLAT PLUG

Valleylab; international models; 3m

unit side Ø 4.0mm; Umax: 0.8 KVp

Instrument-side



Unit-side



SKU: 287-215-30

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Cables](#)
TYPE [Bipolar Cables](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [HF-Cables](#)
TYPE [Bipolar Cables](#)