


## BIPOLAR CABLE FOR THE CONNECTION

of BOWA ARC 350&400 units

4.5m



**SKU:** 287-230-30BR

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Bipolar Cables</a>