


BIPOLAR CABLE FOR THE CONNECTION

of BOWA ARC 350&400 units

4.5m



SKU: 287-230-30BR

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Bipolar Cables