

## BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 1.0mm

15.0cm (6"); flat plug



**SKU:** 022-103-15

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Bipolar Coagulation Forceps](#)  
**Medicon** 88.17.08  
**SPECIFICATIONS** [angled](#), [blunt](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Bipolar Coagulation Forceps](#)  
**Medicon** 88.17.08  
**SPECIFICATIONS** [angled](#), [blunt](#)