

BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 1.0mm

18.0cm (7"); flat plug

SKU: 022-103-18

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
SPECIFICATIONS [angled](#), [blunt](#)