

BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 1.0mm

24.0cm (9.6"); flat plug

SKU: 022-103-24

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bipolar Coagulation Forceps](#)

SPECIFICATIONS [angled](#), [blunt](#)