

## BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 1.0mm

24.0cm (9.6"); flat plug

**SKU:** 022-103-24

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [blunt](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Bipolar Coagulation Forceps](#)

**SPECIFICATIONS** [angled](#), [blunt](#)