

BIPOLAR COAGULATION FORCEPS

straight; blunt; Ø 2.0mm

20.0cm (8"); flat plug



SKU: 022-104-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
Medicon 88.20.16
SPECIFICATIONS [blunt](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Bipolar Coagulation Forceps](#)
Medicon 88.20.16
SPECIFICATIONS [blunt](#), [straight](#)