

## BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 2.0mm

15.0cm (6"); flat plug



**SKU:** 022-105-15

**Categories:** [Forceps](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.17.18
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">blunt</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.17.18
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">blunt</a>