

BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 2.0mm

15.0cm (6"); flat plug

SKU: 022-105-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.18
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.18
SPECIFICATIONS	angled , blunt