

BIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 2.0mm

20.0cm (8"); flat plug

SKU: 022-105-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.20.18
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.20.18
SPECIFICATIONS	angled , blunt