

BIPOLAR COAGULATION FORCEPS

angled on shanks; pointed; Ø 0.3mm

17.0cm (6.8"); flat plug

SKU: 022-107-03

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.25
SPECIFICATIONS	curved , sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.17.25
SPECIFICATIONS	curved , sharp