

## BIPOLAR COAGULATION FORCEPS

angled on shanks; pointed; Ø 0.3mm

17.0cm (6.8"); flat plug

**SKU:** 022-107-03

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.17.25
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">sharp</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.17.25
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">sharp</a>