

BIPOLAR COAGULATION FORCEPS

angled on shanks; blunt; Ø 1.0mm

17.0cm (6.8"); flat plug

SKU: 022-107-10

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.17.26 |
| SPECIFICATIONS | blunt , curved |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.17.26 |
| SPECIFICATIONS | blunt , curved |