

## BIPOLAR COAGULATION FORCEPS

angled on shanks; blunt; Ø 1.0mm

17.0cm (6.8"); flat plug

**SKU:** 022-107-20

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Bipolar Coagulation Forceps](#)  
**SPECIFICATIONS** [blunt](#), [curved](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Bipolar Coagulation Forceps](#)  
**SPECIFICATIONS** [blunt](#), [curved](#)